

# VITA ENAMIC® Polishing Set clinical/technical

For simple pre- and high-gloss polishing of VITA ENAMIC hybrid ceramic



VITA shade taking

VITA shade communication

VITA shade reproduction

VITA shade control

Date of issue: 09.12

VITA shade, VITA made.

**VITA**

## Product description

The VITA ENAMIC polishing sets were developed for reliable, efficient and material-specific surface treatment of hybrid ceramic restorations in dental practices and laboratories. The sets include various polishing instruments for pre- and

high-gloss polishing. These instruments are suitable for careful and gentle polishing of occlusal surfaces, cusps, fissures and contact points of the restoration and produce surfaces with exceptional gloss.

## Product advantages

### Excellent final results

High Gloss and plaque-resistant surfaces can be produced with these instruments. Precise concentricity, matched grit sizes and the individual geometries of the instruments guarantee results with superior precision.

### Simple and safe handling

The instruments provide superior material removal performance. Good handling and the use without polishing paste enable simple and fast processing.

### Gentle and careful processing

These instruments were developed specifically for the VITA ENAMIC to ensure gentle and careful polishing, reducing the risk of possible microcrack formation.

## Product data

VITA ENAMIC Polishing Set clinical: includes a total of 8 instruments for the contra-angle; 4 instruments for pre-polishing and 4 instruments for high-gloss polishing.

VITA ENAMIC Polishing Set technical: includes a total of 8 instruments for the handpiece; 4 instruments for pre-polishing and 4 instruments for high-gloss polishing.

## VITA ENAMIC® sets

Designation	Description	Content / pieces	Prod. No.
VITA ENAMIC Polishing Set clinical	2-stage polishing system with pre- and high-gloss polishers for the contra-angle	8	EENPSETC
VITA ENAMIC Polishing Set technical	2-stage polishing system with pre- and high-gloss polishers for the handpiece	8	EENPSETT

## VITA ENAMIC® Polishing Set clinical – refill packs

Designation	Description	Content / pieces	Prod. No.
<b>Pre-polishing (pink)</b>	SiC instruments (silicon carbide)		
Polishing instrument, small tip, VI-ES3m		6	ERWES3M6
Polishing instrument, tip, VI-ES4m		6	ERWES4M6
Polishing instrument, lens, VI-EL10m		6	ERWEL10M6
Polishing instrument, cup, VI-EC7m		6	ERWEC7M6
<b>High-gloss polishing (grey)</b>	Diamond instruments		
Polishing instrument, small tip, VI-ES3f		6	ERWES3F6
Polishing instrument, tip, VI-ES4f		6	ERWES4F6
Polishing instrument, lens, VI-EL10f		6	ERWEL10F6
Polishing instrument, cup, VI-EC6f		6	ERWEC6F6

## VITA ENAMIC® Polishing Set technical – refill packs

Designation	Description	Content / pieces	Prod. No.
<b>Pre-polishing (pink)</b>	SiC instruments (silicon carbide)		
Polishing instrument, tip, VI-ES5m		6	ERHES5M6
Polishing instrument, torpedo, VI-ET6m		6	ERHET6M6
Polishing instrument, wheel, VI-ER15m		6	ERHER15M6
Polishing instrument, lens, VI-EL16m		6	ERHEL16M6
<b>High-gloss polishing (grey)</b>	Diamond instruments		
Polishing instrument, tip, VI-ES5f		6	ERHES5F6
Polishing instrument, torpedo, VI-ET6f		6	ERHET6F6
Polishing instrument, wheel, VI-ER15f		6	ERHER15F6
Polishing instrument, lens, VI-EL16f		6	ERHEL16F6

More information about VITA ENAMIC is available at [www.vita-enamic.de](http://www.vita-enamic.de) / [www.vita-enamic.com](http://www.vita-enamic.com)

**VITA**