

# System Compatibility



VITA IMPLANT SOLUTIONS blanks with a specific holder system are available for the following CAD/CAM systems:

- CEREC with software version 4.4 (Sirona Dental GmbH)
- inLab with software 15.0 (Sirona Dental GmbH)

VITA IMPLANT SOLUTIONS are compatible\* with the following implant system manufacturers via the integrated interface for the adhesive/titanium base (TiBase, Sirona Dental GmbH):



- Nobel Biocare
- Straumann
- Astra Tech
- Friadent
- Zimmer
- Medentika Implant
- CAMLOG
- Biomet 3i


\* A comprehensive list of implant systems compatible with TiBase is available at [sirona.com](http://sirona.com).

## Geometries

An integrated S or L interface is available for each geometry size.




**IS-14:** 18 x 14 x 12  
(for abutment posts)



**IS-16:** 18 x 16 x 18  
(for abutment crowns)

## Shades



## Ordering Information

Order No.	Patterson Order No.	Translucency	Shade	Qty	Size/Type (L or S)
EC4EM4069485	70492785	HT	1M1	5	16mm L
EC4EM4079485	70492793	HT	1M2	5	16mm L
EC4EM4129485	70492801	HT	2M2	5	16mm L
EC4EM4209485	70492819	HT	3M2	5	16mm L
EC4EM4289485	70492827	HT	4M2	5	16mm L
EC4EM3069725	70492835	T	1M1	5	14mm L
EC4EM3079725	70492843	T	1M2	5	14mm L
EC4EM3129725	70492850	T	2M2	5	14mm L
EC4EM3209725	70492868	T	3M2	5	14mm L
EC4EM3289725	70492876	T	4M2	5	14mm L
EC4EM4069365	70492918	HT	1M1	5	16mm S
EC4EM4079365	70492926	HT	1M2	5	16mm S
EC4EM4129365	70492934	HT	2M2	5	16mm S
EC4EM4209365	70492942	HT	3M2	5	16mm S
EC4EM4289365	70492959	HT	4M2	5	16mm S
EC4EM3069605	70492967	T	1M1	5	14mm S
EC4EM3079605	70492975	T	1M2	5	14mm S
EC4EM3129605	70492983	T	2M2	5	14mm S
EC4EM3209605	70492991	T	3M2	5	14mm S
EC4EM3289605	70493007	T	4M2	5	14mm S

Polishing Sets	Description	Order No.
VITA ENAMIC Polishing Set – clinical	8 pre and high-gloss polishers for contra-angle	EC4EM3209605
VITA ENAMIC Polishing Set – technical	8 pre and high-gloss polishers for contra-angle	EC4EM3289605